


S5[®]

READY TO MEET ANY CHALLENGE



Cardiac surgery solutions



INDEX

	Page
CONCEPT	4 – 7
MAST ROLLER PUMPS	8 – 9
ROLLER PUMPS	10 – 11
CP5 CENTRIFUGAL PUMP	12 – 13
TUBING CLAMPS	14 – 15
SYSTEM PANEL	16 – 17
TIMER	18
PRESSURE CONTROL	19
TEMPERATURE MONITOR	20
LEVEL CONTROL	21
BUBBLE DETECTOR	22
CARDIOPLEGIA CONTROL	23
B-CARE ₅	24
ELECTRONIC GAS BLENDER	25
VENOUS LINE CLAMP	26
EVO - ELECTRICAL VENOUS OCCLUDER	27
ACCESSORIES	28 – 31
HEATER-COOLER 3T	32 – 33
CONNECT	34 – 35
TECHNICAL SPECIFICATIONS	36 – 39

S5™

READY TO MEET ANY CHALLENGE

THE CONCEPT

The S5 heart-lung-machine has a modular design.

The individual components, modules and accessories can be easily replaced, the overall system can be modified and extended, and the modular structure saves valuable time during routine maintenance.



product
design
award

2006 ■



S5 mast roller pump 85 with “C”-shaped holder

This version of the S5 mast roller pump system was specially designed for operators who infrequently perform paediatric and infant/neonate perfusion. The pumps can be placed closer to the patient

leading to a reduction in tubing length and priming volume – a positive effect that is important for paediatric and infant/neonate perfusion, but also for other applications with low flow rates.



CP5 centrifugal pump system for S5

Versatility: the centrifugal pump can be mounted quickly onto a mast and integrated into the available system saving space.

The CP5 uses the monitoring and control units that are already part of the system.



Mast system extension with S5 mast roller pump 150 / two S5 mast roller pumps 85

The S5 roller pumps are mounted on two masts that are connected to an adjustable swivel arm. The control units of the pump heads are mounted separately on the mast

system of the console. Convenient fast clamp connectors allow the mast roller pumps to be mounted/removed quickly.

S5™

CONSOLE

S5 console

The console housing accommodates and protects the entire electronics of the E/P pack including the central power supply and the uninterruptible emergency power supply. The four castors of the console can be locked separately.

S5 pump table

This stainless steel pump table is screwed to the console. Stainless steel pins on the pump table are used for mounting and securing the pump housing. Pump tables for 3, 4 or 5 pumps can be supplied.

S5 standard mast system

The standard mast system is fixed to the console and includes:

Two fixed telescope masts and an adjustable vertical mast with an infusion rack. The masts can be used for mounting the S5 system panel and additional accessories and disposables. The height of all masts can be adjusted.

Two height-adjustable push bars on the left and right side of the console can be used when transporting the S5. They can also be used for mounting accessories.

A horizontal mast stabilises the mast system.



Product designation	Part number		
Consoles with E/P pack and standard mast system	3-position	4-position	5-position
S5 console	48-30-00	48-40-00	48-50-00
Mast systems	Size 3	Size 4	Size 5
Telescope mast with infusion rack	all sizes: 48-30-50		
Telescope mast, movable with infusion rack	all sizes: 48-30-51		
Push bar (horizontal)	all sizes: 48-30-57		
Push bar mast (vertical)	all sizes: 48-30-67		
"C"-shaped holder	all sizes: 50-70-57		
Horizontal mast	48-30-77	48-30-78	48-30-79
Crossbar for movable mast (horizontal)	48-30-81	48-30-82	48-30-83

MAST SYSTEM EXTENSION



S5 Mast System Extension

The S5 mast system extension can accommodate up to 3 mast roller pumps. The separate control units are mounted away from the pump heads on the console mast system. The mast system extension can be mounted on the left or on the right side of the console and folded prior to transport. A supporting castor is used to stabilise the mast system. The mast system extension can be used to place the pumps, oxygenator and tubing set right beside the patient.

The mast system extension

- is for paediatric and infant/neonate perfusion
- can be operated right beside the patient
- can be swivelled and is user friendly



Product designation	Part number
S5 mast system extension	50-45-00
consisting of:	
Swivel telescope mast with infusion rack and castor	50-45-05
2 swivel arms	50-45-10
Vertical mast (including 2 horizontal masts)	50-45-15
Transport locking arm	50-45-20

S5™

MAST ROLLER PUMPS



MAST ROLLER PUMPS

Product designation	Part number
Mast roller pump system 85 consisting of:	50-80-70
1x mast roller pump 85 (connection cable supplied)	10-88-60
Control panel for mast roller pump 150/85 (connection cable supplied)	28-95-80

Product designation	Part number
Mast roller pump system 150 consisting of:	50-80-00
Mast roller pump 150 (connection cable supplied)	10-88-00
Control panel for mast roller pump 150/85 (connection cable supplied)	28-95-80

Product designation	Part number
Mast roller pump system 85 consisting of:	50-80-60
2x mast roller pumps 85 (connection cable supplied)	10-88-60
Control panel for mast roller pump 85 (connection cable supplied)	28-95-85

Product designation	Part number
Mast roller pump system with 2 MRP 85 consisting of:	50-80-62
2x mast roller pumps 85 (connection cable supplied)	10-88-60
Control panel for 2 mast roller pumps 85 (connection cable supplied)	28-95-85

Double holder (fixed)

Two mast roller pumps, with 85 and 150 mm diameters respectively, are available for specialised applications. The small double pump, which is controlled from a single control panel, has some special features.

The holders that were specially developed for the mast roller pump are equipped with practical fast clamp connectors for quick and easy mounting/removal. They can be quickly and easily mounted to/removed from the mast. The corrosion-proof stainless steel pump housing is easily cleaned.

MRP 85

- user friendly
- whenever a robust design and trouble-free upgrading are required

ROLLER PUMPS

S5 roller pump

The S5 roller pump with a 150 mm diameter is used primarily for the arterial blood flow during cardiotomy suction and venting.

S5 double roller pump

The S5 double roller pump combines two roller pumps with a diameter of 85 mm in a single housing. Both roller pumps can be operated and controlled independently. The double roller pump is especially suitable for paediatric/neonate perfusion, as a cardioplegia pump or as a cardiotomy suction and venting pump.

The roller pumps are the most important elements of the S5 perfusion system. Every pump has an independent control system and its own pump control panel that is operated using a high-contrast colour touch screen. The pumps can be individually configured, i.e. the monitoring functions can be individually assigned to each pump and displayed on the touch screen.

An innovative, hard-wearing setting knob (incremental shaft encoder) is used to adjust the set speed electronically. A maximum of five roller pumps can be placed on the console pump table and connected to the E/P pack. The optimised pump head geometry (horse-shoe shaped) keeps pressure peaks during operation to a minimum.



ROLLER PUMPS

Product advantages

- For perfect positioning relative to the surgical site, the pump heads can be rotated 180° (RP 150) or 240° (RP 85). Rotation is carried out in increments of 15° (self-locking).
- The pump head design facilitates quick and easy insertion of the pump tubing.
- The Master-follower mode regulates the flow of the follower pump according to the set flow ratio of the master pump.
- Each pump can be operated in continuous or pulse mode.
- Large digital displays of the speed and/or the flow rate provide a quick overview of the current pump speed.
- Quick orientation: the monitoring functions and current pump status are displayed on the touch screen when an alarm is triggered.
- All pump parameters can be entered using the clearly-structured menu navigation on the pump touch screen.
- If two Override keys are pressed simultaneously, the monitoring functions assigned to a pump are overridden.
- Every pump has a stop link function. This function can be connected to the arterial pump. If the arterial pump stops, the stop link pump also stops.



Product designation	Part number
S5 roller pump 150 / S5 RP 150	10-80-00
S5 double roller pump 85 / S5 DRP 85	10-85-00

CP5 CENTRIFUGAL PUMP

This established centrifugal pump system provides the S5 operator with flexibility. The mast assembly saves space, no additional space is needed on the console.

The ERC closes the arterial line in a fraction of a second if an alarm is triggered (e.g. retrograde flow, level or bubble alarm), keeping the danger of air delivery to a minimum.

The overall system is delivered complete with drive unit, control panel, connection cable, sensor module for flow measurement, flow sensor and emergency system.





Product designation	Part number
---------------------	-------------

CP5 (only for S5 C5) consisting of:	60-00-60
---	----------


CP5 drive unit (with cable)	60-01-04
-----------------------------	----------

Pump control panel (with holder)	60-02-60
----------------------------------	----------

Emergency system (complete)	60-01-35
-----------------------------	----------

Flow sensor (3/8")	96-414-140
--------------------	------------

Flow measurement sensor module	25-60-70
--------------------------------	----------



Product designation	Part number
---------------------	-------------

Electrical remote-controlled tubing clamp consisting of: tubing clamp, 3-joint mast holder with fast clamp connector, 500 mm	60-05-60
--	----------

Electrical remote-controlled tubing clamp consisting of: tubing clamp, 3-joint mast holder with fast clamp connector, 620 mm	60-05-65
--	----------

TUBING CLAMPS

Tubing clamp inserts are available for all tubing sizes. They are inserted into the tubing clamp block of the pump heads.

Special cardioplegia tubing clamp inserts allow two tubes with different diameters to be safely and simultaneously connected to the roller pump heads. They are available for flow ratios from 1:1 to 8:1.




Product designation				Part number
Tubing clamp block RP 150 (included in RP 150, 10-80-00)				10-81-35
Tubing clamp inserts for tubing clamp block RP 150				
	1/4" x 1/16"	red		10-64-15
	1/4" x 3/32"	yellow		10-64-25
	5/16" x 1/16"			
	3/8" x 1/16"	black		10-64-40
	5/16" x 3/32"			
	3/8" x 3/32"	blue		10-64-50
	1/2" x 3/32"	green		10-64-65
	1/8" x 1/16"	violet		10-64-05
	3/16" x 1/16"	mint-green		10-64-10
	1/2" x 1/16"	grey		10-64-55
Tubing clamp inserts for cardioplegia RP 150				
	3/16" x 1/16"	light grey	1:1	10-64-70
	3/16" x 1/16"			
	1/4" x 1/16"	light brown	1:1	10-64-71
	1/4" x 1/16"			
	3/16" x 1/16"	white	2:1	10-64-72
	1/8" x 1/16"			
	1/4" x 1/16"	light blue	2:1	10-64-74
	3/16" x 1/16"			
	1/4" x 1/16"	turquoise	4:1	10-64-76
	1/8" x 1/16"			
	17/64" x 1/16"	brown	8:1	10-64-78
	3/32" x 1/16"			

TUBING CLAMPS

Variolock tubing clamp module
The Variolock tubing clamp module was developed for the heads of roller pumps. A large range of tubing clamp

inserts for single and double (e.g. for cardioplegia delivery) tubing configuration are available.



Product designation			Part number
Variolock tubing clamp module RP 150 (optional)			10-81-30
Tubing clamp inserts for Variolock			
included (in 10-81-30)	1/4" x 1/16"	red (small)	10-61-73
	1/4" x 3/32"		
	5/16" x 1/16"		
	3/8" x 1/16"		
	5/16" x 3/32"		
	3/8" x 3/32"		
3/8" x 3/32"	blue (large)	10-61-72	
	1/2" x 1/16"		
	1/2" x 3/32"		
Tubing clamp inserts for Variolock for cardioplegia			
optional	1/4" x 1/16"	1:1	10-61-91
	1/4" x 1/16"		
	1/4" x 1/16"	2:1	10-61-92
	3/16" x 1/16"		
	1/4" x 1/16"	4:1	10-61-93
	1/8" x 1/16"		

Tubing clamp block DRP 85


Tubing clamp blocks that can take tubing sizes up to 5/16" x 1/16" have been developed for the small roller pump.



Product designation			Part number
Tubing clamp block DRP 85 (included in DRP 85, 10-85-00)			
included	1/4" x 1/16"	red	10-86-56
	1/4" x 3/32"	yellow	10-86-57
	5/16" x 3/32"	black	10-86-58
optional	1/8" x 1/16"	violet	10-86-55
	3/16" x 1/16"	mint-green	10-86-59

S5™

SYSTEM PANEL



S5 system panel

The S5 system panel contains the display and control modules for all of the monitoring, control and measuring devices and is, alongside the pump control panel, another interface between the operator and the S5 System. The system panel can be mounted on the left or right mast of the standard mast system as required.

The holder with ball joint can be swivelled into any desired position.

System panels with 3 to 6 slots accommodate the control and monitoring modules. All display and control modules can be replaced during operation if a fault occurs. The data is displayed unchanged on the replacement module.

SYSTEM PANEL

**S5 blank module**

If not all of the slots are taken up by display and control modules and adding additional monitoring parameters should remain an option, the inexpensive blank module can be used as a placeholder.

**S5 display and control module**

All display and control modules are physically identical but each one is controlled by its own separate microprocessor. The high-contrast TFT display has a restrained colour scheme and a clear display structure. The display layout is determined by the chosen control and monitoring function. All pump control and monitoring function settings - with a few exceptions - are entered using the control module touch screens. A user-friendly menu guides the operator through all setting parameters.


Product designation	Part number
S5 system components	
3-position S5 system panel for 3 display and control modules	28-95-03
4-position S5 system panel for 4 display and control modules	28-95-00
5-position S5 system panel for 5 display and control modules	28-95-01
6-position S5 system panel for 6 display and control modules	28-95-04
1 display and control module	28-95-10
Blank module	28-95-30

TIMER


Three timers that work independently of each other can measure the duration of three simultaneous separate processes, for example the complete bypass time or the aortic cross-clamping time.

- The timers can be started and stopped individually.
- It is possible to carry out cumulative measurements with each timer.
- The measurement range for each timer is 999 min 59 sec.

A fourth timer can (depending on the setting) count upwards or downwards for a maximum of 10 hours or 600 minutes (optional, available on request).



PRESSURE CONTROL




The double pressure sensor module is used for:

Measuring and displaying the pressure in the extracorporeal circuit. The display range extends from -200 mmHg to +800 mmHg. The values can be displayed in either mmHg or kPa.

Limiting the pressure by stopping the pump (monitoring mode) when the preset pressure (stop limit) has been reached.

Controlling perfusion with constant adjustable pressure (set value) through automatic variation of the pump speed (control mode).

The double pressure sensor module allows the operator to control two pumps independently of one another. The control and display-module of the system panel is used for setting parameters and adjusting the zero-point. The pressure transducer and transducer holder are required accessories.





Product designation	Part number
Sensor module 2-channel pressure monitor	22-20-20
Sensors and accessories (optional)	Part number
Medex transducer (MX 960)	45-04-03
Cable for Medex transducer	45-04-15
Holder for one Medex transducer	45-04-16
Holder for 2 Medex transducers	45-04-17




TEMPERATURE MONITOR

- The temperature monitor allows the simultaneous measurement and display of up to four temperatures. One temperature sensor module channel is reserved for cardioplegia control.
- An upper and a lower temperature limit can be set. The values are set on the control and display module on the system panel.
- An outstanding feature of the measuring channels is their measurement precision.
- If the temperature limit is reached, visual and acoustic alarms are triggered.
- Temperature probes (listed under "Accessories") are required to operate the monitor.



Product designation	Part number
Sensor module 4-channel temperature monitor	20-30-20
Accessories (optional)	Part number
Temperature probes for direct measurement in the oxygenator (optional)	
For SORIN GROUP oxygenators	45-03-10
For Dideco oxygenators	45-03-11
For MEDTRONIC oxygenators	45-03-30

LEVEL CONTROL




- The level monitor controls the blood level in the oxygenator/reservoir.
- The level sensor is mounted on special disposable, adhesive mounting strips which are attached to the reservoir or oxygenator. The sensor can be used with most hardshell oxygenators and reservoirs.
- The sensor is based on the effect the liquid level has on high-frequency, electromagnetic waves in an electric oscillating circuit. Using this method, interference caused by electromagnetic or external light sources and residual liquid on the oxygenator's wall is avoided.

Two operating modes can be used for level monitoring and control:

- In monitoring mode, the sensor is attached to the minimum filling level. If the current blood level drops below the set minimum level, a visual and acoustic alarm is triggered and the pump stops. The level icon switches to red. As soon as the blood level rises above the minimum level, the alarm is cleared automatically and the pump restarts. The level icon turns green.
- In control mode, the sensor is attached to the set filling level. By automatically controlling the pump's rotational speed, the set blood level in the reservoir can be kept constant. The display is yellow.

Product designation	Part number
Sensor module level control	23-40-00
Sensors and accessories (included)	Part number
Level sensor	23-27-40
Adhesive sensor holder(100 units)	23-27-41
Accessories (optional)	Part number
Level sensor holder	75-521-548



BUBBLE DETECTOR

This monitoring device detects air bubbles and micro-bubbles in the extracorporeal circuit.

This function requires a sensor module and the corresponding sensor. The scope of delivery includes a bubble detector that can be used for 3/8" tubing. Alternatively, 1/2" and 1/4" sensors are available.

- The 3-joint mast holder with a fast clamp connector allows the bubble sensor to be positioned at the tubing system.
- If air bubbles are detected, a visual and acoustic alarm is triggered and the connected pump stops.
- You can set three different bubble detection alarm thresholds (4 mm, 5 mm and 6.5 mm Ø).
- The micro-bubble detection function can also be activated.



Product designation

Part number

Sensor module bubble detector consists of:

One bubble sensor 3/8" (9.56 mm) (23-07-50)
and one 3-joint mast holder, 420 mm (23-26-96)

23-45-00

Sensor module bubble detector consists of:

One bubble sensor 1/4" (6.35 mm) (23-07-40)
and one 3-joint mast holder, 420 mm (23-26-96)

23-45-01

Sensor module bubble detector consists of:

One bubble sensor 1/2" (12.7 mm) (23-07-45)
and one 3-joint mast holder, 420 mm (23-26-96)

23-45-02

Sensor module bubble detector consists of:

One bubble sensor 3/8" (9.56 mm) (23-07-50)
and one 3-joint mast holder, 620 mm (23-26-91)

23-45-10

Sensor module bubble detector consists of:

One bubble sensor 1/4" (6.35 mm) (23-07-40)
and one 3-joint mast holder, 620 mm (23-26-91)

23-45-11

Sensor module bubble detector consists of:


One bubble sensor 1/2" (12.7 mm) (23-07-45)
and one 3-joint mast holder, 620 mm (23-26-91)

23-45-12

Ultrasound gel (required for 3/8" or 1/2") 250 ml

96-06-10

CARDIOPLEGIA CONTROL



This unit can be used with a RP 150 or a DRP 85 to deliver cardioplegic solutions or blood cardioplegia during an operation. The operator can choose between two operating modes that can be selected in the menu of the control and display module.

Manual operation

The operator can start and stop the pump. The dose volume to be delivered counts up on the volume display (beginning at 0).

Automatic operation

In this operational mode an exact preset dose is delivered. In this case the volume display starts at the preset dose and counts down to 0. Then the pump stops.



- The integrated timer automatically starts during a pump stop regardless of operational mode and records the ischaemia time.
- The (total) volume delivered since the start of cardioplegia is accumulated and displayed.
- Two roller pumps are particularly suited to blood cardioplegia delivery.
- The cardioplegia sensor module has its own connectors for a bubble sensor and a pressure transducer.
- As soon as the sensor detects bubbles, the cardioplegia pump stops automatically and the cardioplegia delivery is interrupted. At the same time, the visual and acoustic alarms are triggered.
- When the preset pressure (stop limit) has been exceeded, the cardioplegia pump stops and cardioplegia delivery is interrupted. At the same time, the visual and acoustic alarms are triggered (monitoring mode). The control mode can also be set.

Product designation	Part number
Sensor module cardioplegia control	27-80-20
Sensors and accessories (optional)	Part number
Medex transducer (MX 960)	45-04-03
Cable for Medex transducer	45-04-15
Holder for one Medex transducer	45-04-16
Holder for 2 Medex transducers	45-04-17
Bubble sensor	
1/4" (6.35 mm)	23-07-40
1/2" (12.7 mm)	23-07-45
3/8" (9.56 mm)	23-07-50
Ultrasound gel (only required for 23-07-45 or 23-07-50), 250 ml bottle	96-06-10
3-joint mast holder with fast clamp connectors for 2 sensors, 620 mm	23-26-91
3-joint mast holder with fast clamp connectors for 2 sensors, 420 mm	23-26-96

B-CARE₅

The In-Line Monitoring System for 3 vital parameters

B-Care₅ is the first in-line monitor fully integrated into a heart-lung-machine.


Vital parameters, i.e., venous saturation, haematocrit and venous temperature are measured and shown on a display and control module (DCM) of the S5 system panel. This new S5 sensor module is plugged into the E/P pack of the S5 perfusion system, eliminating the need of an external blood monitor. Moreover, all measured blood values are automatically transferred to the Data Management System via the S5's internal data communication system.

Venous Connector

The venous connector is available either separately packaged and sterile or pre-connected in a customised tubing system.


Venous Sensor

B-Care₅ performs an automatic self-test upon power-up of the S5 (no intervention required by the user). The venous sensor features an optical sensor to measure signals determining the venous saturation (Sat %) and the haematocrit (Hct %). An integrated thermistor reads the actual venous blood temperature.



Product designation	Part number
B-Care ₅	25-60-00
Sensors and accessories (optional)	Part number
B-Care ₅ sensor module (venous)	25-60-20
Venous sensor (complete)	97-231-059
Holder for measuring head (reference element)	25-60-30
Connectors	
1/2" venous measuring connector 10 units (sterile)	05171
3/8" venous measuring connector 10 units (sterile)	05172
Kit for 1/4" venous connector	080099

ELECTRONIC GAS BLENDER



The electronic gas blender allows to set, monitor and display the gas flows required for extracorporeal circulation. The preset values (i.e. the total flow including Air + O₂, FiO₂ and CO₂) can be set independently and are displayed on both the gas blender and the display and control panel.

The actual values and the set values are continuously compared. Additionally, the actual value is measured by 2 independent sensors and an alarm is triggered if a deviation between the 2 values is detected.

The remote control on the display and control module can be used to change the set values for Air + O₂, FiO₂ and CO₂ from gas flow to blood flow. The operator is made aware of the actual values exceeding or dropping below the set values by acoustic and visual signals.

The electronic gas blender is available in three different versions:

- Electronic gas blender (10 l/min) for adult perfusion
- Electronic gas blender (5 l/min) for paediatric perfusion
- Electronic gas blender (2 l/min) for infant/neonate perfusion

Product designation	Part number
Electronic gas blender (10 l/min)	25-28-67
Electronic gas blender (5 l/min)	25-28-68
Electronic gas blender (2 l/min)	25-28-69
Standard holder (straight, with fast clamp connector)	55-91-50
Accessories (optional)	Part number
Holder (U-shaped with plate, only compatible with 4- and 5-position consoles)	25-40-70

Scope of delivery:

- Electronic gas blender
- Connection cable for E/P pack (length: 2 m)
- 4 connectors (Air, O₂, CO₂, total gas flow)

VENOUS LINE CLAMP

Venous Line Clamps with Mechanical Remote Control

The venous line clamp has a mechanical remote control.

It has a lightweight design. If the clamp head is fixed to another part of the venous tubing system, it might well be necessary to support the clamp head with a joint holder.

The 1 m Bowden cable connects the clamp head with the control unit and transfers settings entered on the control unit immediately to the clamping lever in the clamp head. The control unit is mounted on one of the S5's push bars (right or left – depending on the system arrangement and/or ease of use). The coarse and fine setting knobs are used for setting the tubing diameter and regulating the venous return flow quickly.




An overview of the remote-controlled venous line clamp's advantages:

- The separate control unit and line clamp facilitates a flexible assembly.
- Unlike an electronic line clamp, this clamp does not have to be calibrated and is always ready for use.
- The structural design of the line clamp means that the clamped line can be visually monitored at all times.

Product designation	Part number
Venous line clamp remote control	12-40-00
3-joint mast holder for venous line clamp	Part number
With fast clamp connector (586 mm)	12-30-90
With fast clamp connector (386 mm)	12-30-95
Mast adapter for mounting the control unit on the vertical mast	12-05-80
Tubing Inserts in sets with 4 pcs. each (incl.)	Part number
∅ 1/4" x 1/16" red	10-07-20
∅ 3/8" x 3/32" blue	10-07-23
∅ 1/2" x 3/32" green	10-07-25
Tubing Inserts in sets with 4 pcs. each (opt.)	Part number
∅ 1/4" x 3/32" yellow	10-07-21
∅ 3/8" x 1/16" black	10-07-22
∅ 1/2" x 1/16" grey	10-07-24
∅ 5/8" x 3/32" brown	10-07-26
∅ 1/8" x 1/16" violet	10-07-27
∅ 3/16" x 1/16" turquoise	10-07-28

EVO – ELECTRICAL VENOUS OCCLUDER



The clamp closes automatically when the stop link function to the arterial pump is activated, if the latter has been stopped by monitoring functions in case of an alarm or is stopped manually. When the arterial pump starts up, the EVO opens to the most recently specified set value. An override of the stop link function is possible at any time, directly at the EVO operating unit.

Use the relevant keys or the setting knob on the EVO operating unit to open and close. When the setting knob is turned, there is audible clicking and locking into place. Different ranges can be selected for a fine adjustment.


- Ergonomic operation when initiating and ending ECC
- The set value can be preset when the occluder is closed and the stop link function is activated
- Audible clicking and locking into place make you doubly aware of setting knob adjustments
- A choice of fine adjustment from < 40% in 10% steps
- The stop link delay for a level alarm is adjustable from 0 to 60 seconds

Product designation	Part number
Electrical venous occluder consists of:	12-80-00
Occluder (with mast holder)	12-80-10
Control unit	28-95-70

ACCESSORIES

Colour coding set

To ensure the control panel is safely assigned to the mast roller pump, a colour coding set is available.



Cable holders

Cable holders are used for the secure and correct routing of cables and tubing on the mast system. The holders can be installed anywhere on the mast system. Available in sets of 6.



Perfusionist's chair

This comfortable chair helps the operator maintain a healthy posture. The ergonomic swivel design is mounted on 5 self-locking castors. The back and seat height are adjustable. It is suitable for the operating theatre and easy to clean.



S3 to S5 adapter cable

The adapter cable is used for connecting S3 accessories (SCP, ERC, gas blender) to the S5 System.



Product designation	Part number
Colour coding set	50-80-99
Cable holders / 6 clips Ø 33 mm mast diameter	45-09-10
Cable holders / 6 clips Ø 25 mm mast diameter	45-09-11
Perfusionist's chair	41-02-98
S3 to S5 adapter cable	45-12-00

ACCESSORIES

**LED console lamp**

A second version of the console lamp uses LEDs. Some major advantages of LEDs are:

- Powerful luminosity
- Longevity
- Energy efficiency
- Minimum heating of lamp

**Writing desk DIN A4**

This writing desk is mounted on the console using a stable arm. The position and tilt can be set individually by the operator.

Cover for fast clamp connectors

These convenient covers safeguard the fast clamp connectors against unintentional opening and protect against spilled liquid.



Product designation	Part number
LED console lamp	35-05-80
Writing desk DIN A4	48-04-00
Cover for fast clamp connectors (set of 6 pieces)	43-42-61

ACCESSORIES

S5 shelf

This versatile shelf is available in three sizes (for 3-, 4- and 5-position pump tables). It can be supplied with or without a AC outlet strip.



S5 pump spacer

The S5 roller pumps for 3-, 4- and 5-position pump tables can be raised by 10 cm using the pump spacer.



S5 plexiglas display protector

This cover protects the roller pump displays against falling objects.



S5 shelf for roller pump 150

A stainless steel shelf for storing small parts is available for every roller pump. A longer version for 2 roller pumps is also available. For technical reasons, a shelf is not available for the double roller pump.



Product designation	Part number		
	3-position	4-position	5-position
S5 shelf with AC outlet strip	48-31-19	48-41-19	48-51-19
S5 shelf without AC outlet strip	48-31-10	48-41-10	48-51-10
S5 pump spacer	48-31-20	48-41-20	48-51-20
S5 plexiglas display protector	48-31-30	48-41-30	48-51-30
S5 shelf for roller pump 150	all sizes: 10-84-60		
S5 shelf for 2 roller pumps 150	all sizes: 10-64-84		


ACCESSORIES

**S5 drawer module**

The stainless steel drawer module is used for storing utensils. It has an easy-access drawer on tracks with a stop to prevent it from falling out. There is an additional sliding tray inside the drawer. The drawer module occupies a single pump space on the console.

**S5 ice container**

Infusion bags and bottles (for example cardioplegic solutions) can be cooled and stored in the S5 ice container. It consists of an outside casing and a stainless steel insert. The chilled bottles etc. are always within reach if the ice container is mounted on the right or left push bar of the S5!

**S5 tubing guide holder**

The tubing guide consists of a 6-m length of PVC tubing with connectors and 3 holders for securing to the horizontal crossbar of the movable mast. The holder lifts the tubing clear of the floor so that it is not squashed under the castors of the console. The floor is also easier to clean without tangled tubing and the danger of being tripped.

Product designation**Part number**

S5 drawer module

48-41-70

S5 ice container

16-05-40

S5 tubing guide holder


16-05-60

incl. connectors and a 6-m length of PVC tubing

HEATER-COOLER 3T

The Heater-Cooler 3T facilitates the precise and fast regulation of the patient's blood temperature and the temperature of the cardioplegic solution regardless of the hot/cold water supply in the operating theatre.

The device has 3 separate water tanks and 3 water circuits that can be used simultaneously. Circuits 1 and 2 use an identical preset temperature and are mainly used to control the patient's temperature.



Circuit 3, which has a separate cold water and warm water tank, is specially designed for cooling and heating blood and/or cardioplegic solutions. Cold water and warm water tanks with the relevant preset temperatures are available during operation at all times.

Heating-cooling blankets connected to the heater-cooler provide an additional support for regulating the patient's blood temperature. The device is operated and monitored from its own control panel or, alternatively, from the display and control modules on the system panel (see illustration).

Product advantages:

- The low tank volume makes short pre-cooling/pre-heating times possible.
- Separate cold and warm water tanks allow the operator to switch between warm and cold cardioplegia cardioplegia spontaneously.
- The patient and cardioplegia circuits can be switched off separately when not in use. This improves the other circuit's heating and cooling performance.
- An independent safety system stops the water temperature reaching critical values.
- The pump suction stage ensures that the heat exchanger and tubing are purged.

HEATER-COOLER 3T



Product designation	Part number
Heater-cooler 3T, 230 V	16-02-80
Accessories (included)	Part number
Cable for connection to the S5 System, 6-m	45-12-16
Accessories (optional)	Part number
Heating-cooling blanket for adults (55 x 150 cm)	16-10-50
Heating-cooling blanket for children (55 x 90 cm)	16-10-51
Tubing set for heating-cooling blanket	16-10-55
1 adapter connector set for blanket	16-10-10



CONNECT

CONNECT is an innovative and intuitive Perfusion charting system focused on real time and retrospective calculations and trending tools to assist with data management during and after CPB.

CONNECT offers a combination of intelligently designed software and hardware to support your perfusion and documentation goals. CONNECT ensures data integrity for an efficient, effective and automated electronic charting system.

The CONNECT Manager application manages all case data in one central database and allows access to the information by the Perfusion team at anytime. Manager can be used for printing or exporting case reports as well as for conducting data analysis with the included statistics tool.


The CONNECT Recorder application effortlessly collects customized patient and case information during CPB. Information from the heart lung machine and other devices in the operating room are displayed in a variety of chart and table formats to optimize viewing preferences. Comments, drugs and volumes are easily entered into the case record via quick entry touch of the innovative tag cloud feature.


CONNECT Recorder operates on a touch screen datapad, custom designed as a drip proof, medical grade tablet PC that mounts onto the Sorin heart lung machine. The touch screen interface enables fast and easy case event entries and patient specific documentation via on screen virtual keyboards, pull down lists and customized data tables. The datapad standard configuration includes built in wireless and LAN connections for transfer of case data and customized settings across the hospital network.

Once the case has been completed, the case record may be transferred between Manager and Recorder via USB drive, across the hospital network with a LAN connection or by utilizing the integrated wireless feature. The record can then be printed, exported to a location on the hospital network and retained in the database for statistical analysis, export and reporting at anytime.

Designed with Goal Directed Perfusion in Mind

- GDP Monitor feature for real time calculation of DO₂, VO₂ and VCO₂ values
- Dedicated fluid balance chart for focus on the hemodilution footprint
- Embedded Quality Scores feature for practice improvement initiatives
- Electronic connection to LivaNova ATS devices to integrate ATS information into the Perfusion record
- Clinical statistics and Inventory reporting tool
- Paperless or printed user defined case reports
- Easy customization for tailoring Connect to Perfusion practice preferences
- Customizable charts and tables for real time and retrospective data trending and analysis
- Connect is a registered medical device





Specifications

CONNECT Manager

- Operating system: Microsoft® Windows® XP service pack 3 / Microsoft® Windows® 7 Ultimate 32 bit
- Random access memory (RAM): 1 GB (recommended)
- Screen resolution: 1024 x 768 pixels
- Database used: Microsoft® SQL Server 2008 R2 Express (recommended) and higher
- Free hard disk space - software: 120 MB / Microsoft® SQL Server 2008 R2

Datapad for CONNECT Recorder

- TFT Flat Panel Display with LED retro-illumination: 4:3 aspect ratio, diagonal of active area 15.0", XGA resolution (1024x768 pixels)
- Operating system: Microsoft® Windows® 7 Ultimate 32 bit
- Intel Celeron 2002E 1.5 GHz
- SO-DIMM DDR3L 1600-RAM 4 GByte
- SSD 2.5" 64 GByte
Removable CFAST 2.0 16 GByte
- 1x COM Port (COM1), TTL Level
- 4x USB 2.0 Ports (EHCI)
- 1 Gigabit LAN 10/100/1000 BaseT controller
- 1 DVI Output port
- WLAN: WPEA-121N 802.11a/b/g/n Half Mini PCIe Module, Atheros AR9382, 2T2R

Product designation	Part number
Connect Manager	24-90-45
Connect Recorder	24-90-80

Additional packages for upgrades from DMS are available, please contact your local LivaNova representative for more details.

1.0 DIMENSIONS, WEIGHTS, OPERATING CONDITIONS

TECHNICAL SPECIFICATIONS

1.1 CONSOLE

Height (to the surface of the pump cover)	640 mm
Depth	600 mm

Console	3-position	4-position	5-position
Width (incl. push bars)	745 mm	890 mm	1073 mm
Weight	83.4 kg	86.3 kg	89.5 kg

Operating conditions	
Operating temperature	+ 10 °C ... + 40 °C
Storage temperature	0 °C ... + 40 °C
Relative humidity (operating and storing)	30% ... 75%

1.2 MASTS

Maximum permissible load	
Maximum total load on mast system	45 kg
Maximum load on a mast	20 kg ⁽¹⁾

Mast system extension (optional)	
Maximum load on the telescope mast	40 kg
Maximum load on the vertical mast	11.5 kg

1.3 PUMPS

	Roller pump 150	Double roller pump 85	Mast roller pump 150	Mast roller pump 85	2 mast roller pumps 85
Height	285 mm	257 mm	289 mm	237 mm	237 mm
Width	180 mm	180 mm	178 mm	116 mm	260 mm
Depth	485 mm	485 mm	299 mm	175 mm ⁽³⁾	200 mm ⁽⁵⁾
Weight	15 kg	12 kg	11.9 kg ⁽²⁾	5 kg ⁽²⁾	11 kg ⁽⁴⁾

Pump specifications	Roller pump	Double roller pump
Diameter of pump raceway Ø	150 mm	85 mm
Diameter of occlusion roller Ø	30.5 mm	15 mm

Speed range	0 to 250 rpm (clockwise, counterclockwise)	
Deviation in speed accuracy	±1% of the terminal value 250 rpm plus ±0.5% of set value	
Speed deviation in the event of a fault (Detection of faulty speed from 30 rpm)	during continuous operation: +15% max.; 2 revolutions max. until pump stops	
Direction of rotation	Clockwise/counterclockwise	Clockwise/counterclockwise

Concentricity		
Pump raceway	0.03 mm	0.03 mm
Occlusion symmetry	0.03 mm	0.03 mm
Occlusion rollers	0.015 mm	0.015 mm

⁽¹⁾ max. swivel arm 200 mm; ⁽²⁾ with fast clamp connector; ⁽³⁾ without fast clamp connector;

⁽⁴⁾ with double holder; ⁽⁵⁾ without double holder

TECHNICAL SPECIFICATIONS

Displays	Roller pump	Double roller pump
rpm display range	0 to 250 rpm	0 to 250 rpm
Resolution	1 rpm	1 rpm
l/min display range (flow)		
1/8"	0 to 0.83 l/min	0 to 0.44 l/min
3/16"	0 to 1.79 l/min	0 to 0.93 l/min
1/4"	0 to 3.12 l/min	0 to 1.57 l/min
5/16"	0 to 4.70 l/min	0 to 2.33 l/min
3/8"	0 to 6.50 l/min	0
1/2"	0 to 11.2 l/min	0

Deviation of speed slave pump	max. 1 percentage point of the flow ratio setting
-------------------------------	---

Power supply	Roller pump	Double roller pump
Operating voltage	24 V DC	24 V DC
Power consumption	160 W	160 W

1.4 SYSTEM PANEL

	For 3 display and control modules	For 4 display and control modules	For 5 display and control modules	For 6 display and control modules
Height	475 mm	590 mm	723 mm	475 mm
Width	184 mm	184 mm	184 mm	375 mm
Depth (without mast holder)	94 mm	94 mm	94 mm	94 mm
Weight (without display and control module)	3.9 kg	4.5 kg	5.1 kg	7 kg

	Display and control module	Control module for mast roller pumps
Height	125 mm	260 mm
Width	179 mm	190 mm
Depth 4	8 mm	100 mm
Weight	0.5 kg	3.5 kg ⁽⁶⁾

2. ELECTRICAL SPECIFICATIONS

2.1 ELECTRONICS AND POWER PACK

Input voltages	100 V ~ bis 240 V~; 50 / 60 Hz
Permissible mains voltage fluctuation	± 10%
Maximum power consumption (standard equipment)	1000 W

⁽⁶⁾ with holder

TECHNICAL SPECIFICATIONS

2.2 UPS and Batteries

Operating time of UPS	
At 400 W output power	20 minutes
At 160 W output power	90 minutes
Charging time	12-15 hours

2.3 System Panel

Display and control module / touch screen	
Operating voltage	24 V
Power consumption	45 W
Pixel Failure Class	Conformity with Pixel Failure Class III

2.4 Shelf with AC Outlet

3-/4-/5-position	
Weight – shelf	approx. 6.5 kg
Maximum load – shelf	8 kg
Number of sockets	4
Protection	at 230/240 V: Circuit breaker 2 A at 110/115 V: Circuit breaker 2 A
Load rating	2 A maximum in total
Sum of leakage currents	500 µA max. in total

Level	
Alarm limit (level sensor) for oxygenators/reservoirs made of rigid polycarbonate, wall thickness at sensor position 3 mm max.	Level display of the sensor holder ± 10 mm

Pressure	
Measurement range mmHg	-200 mmHg to +800
Resolution	1 mmHg

Cardioplegia	
Pressure measurement range	-200 mmHg to +800 mmHg
Resolution	1 mmHg

Temperature monitor	
Display range	0 °C to +50 °C

Timer	
Counting range	0 - 999 min 59 sec

Timer (optional)	
Counting range	0 - 10 h (up and down)

TECHNICAL SPECIFICATIONS

2.5 Modules and Sensors

Level sensor module	
Alarm limit (level sensor) for oxygenators/reservoirs made of rigid polycarbonate, wall thickness at sensor position 3 mm max.	Level display of the sensor holder ± 10 mm

Sensor module for bubble detector	
Alarm limit (bubble sensor) at ≥ 15 rpm 1/2" and 3/8"	Air volume: 0.144 cm ³ (\varnothing 6.5 mm) Air volume: 0.065 cm ³ (\varnothing 5.0 mm) Air volume: 0.034 cm ³ (\varnothing 4.0 mm)

Sensor module 2 channel pressure monitor	
Accuracy	± 5 mmHg
Zero point adjustment range	± 100 mmHg
Gain adjustment range (matching)	$\pm 20\%$
Input resistance	100 k Ω
Output voltage to pressure transducer	< 10 V

Cardioplegia sensor module	
Volume control	
Setting range	0 to 2 liter
Accuracy of dosage	$\pm 10\%$, min. ± 20 ml
Pressure monitor	See sensor module 2 channel pressure monitor
Bubble detector	See alarm limit of the bubble sensor

Sensor module 4 channel temperature monitor	
Temperature measurement range	0 °C to +50 °C
Resolution	0.1 °C
Accuracy (without sensors)	0.0 °C - 25.0 °C ± 0.2 °C 25.0 °C - 45.0 °C ± 0.1 °C 45.0 °C - 50.0 °C ± 0.2 °C

LivaNova

Health innovation that matters

www.livanova.com



LivaNova Deutschland Quality System complies with:
EN ISO 13485:2012

CE⁰¹²³ According to Annex II (Full Quality System) of
MDD 93/42/EEC as amended by directive 2007/47/EEC

Manufacturer:
LivaNova Deutschland GmbH
Lindberghstraße 25
D-80939 München Germany
T.: +49(0)89.32301.0
F.: +49(0)89.32301.555

© 2018 LivaNova all rights reserved.