


NEW

LivaNova

Health innovation that matters

# PureFlex<sup>®</sup> Arterial Cannulae



Cannulae that

can

A complete range of adult arterial cannulae crafted to meet diverse clinical needs in conventional cardiac surgery



## New, Curved and Straight Tip Arterial Cannulae

The all new family of curved and straight tip PureFlex Arterial Cannulae have been designed to meet diverse clinical needs with different configurations. These cannulae provide smooth no-jet flow and a low pressure drop\*, while the new tip design grants ease-of-insertion. The upgraded, sleek DEHP-free tubing is highly transparent for greater visibility and its new wire-reinforced, flexible spring design, with reduced pitch, allows no kinking.





## Fits all Adult Patient Sizes

The PureFlex Arterial Cannulae can boast 4 different sizes and 36 different configurations all of which with a 3/8" tubing connection and is designed to fit all adult patients. The transparent, straight and curved tips come with either a suture flange or collar to support all types of insertion procedures, while the more generous clamping area assures easier clamping during priming.



**PureFlex Arterial Cannulae**  
18, 20, 22, 24 Fr sizes available


## Pressure Drop\*



\* All data extracted from relative Instructions For Use

## Three Different Connections; Greater Versatility


Without luer lock · With luer lock for venting · With porous vent plug and no connector



### Available Configurations


#### Bent tip with collar

RAXX-1BB




Bent tip with collar, with connector 3/8" with luer lock

RAXX-1BC



Bent tip with collar, with connector 3/8" without luer lock


RAXX-1BP



Bent tip with collar, with 3/8" connection and vent plug

#### Bent tip with suture flange

RAXX-2BB



Bent tip with suture flange, with connector 3/8" with luer lock

RAXX-2BC



Bent tip with suture flange, with connector 3/8" without luer lock


RAXX-2BP



Bent tip with suture flange, with 3/8" connection and vent plug

#### Straight tip with collar

RAXX-1SB



Straight tip with collar, with connector 3/8" with luer lock

RAXX-1SC



Straight tip with collar, with connector 3/8" without luer lock

RAXX-1SP



Straight tip with collar, with 3/8" connection and vent plug

# LivaNova

Health innovation that matters

Discover more at [cannulae.livanova.com](http://cannulae.livanova.com)

[www.livanova.com](http://www.livanova.com)



CE 0344

## Sorin Group Italia Srl

A wholly-owned subsidiary of LivaNova PLC  
Via Statale 12 Nord, 86 - 41037 Mirandola (MO) Italy  
Tel: +39 0535 29811 - Fax: +39 0535 25229  
[info.cardiacsurgery@livanova.com](mailto:info.cardiacsurgery@livanova.com)

Please always refer to the Instructions For Use (IFU) manual provided with each product for detailed information, warnings, precautions and possible adverse side effects.

© 2018 LivaNova all rights reserved.